

Product datasheet for **MR206133L4V**

Bcat2 (NM_009737) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bcat2 (NM_009737) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bcat2
Synonyms:	Bcat(m); Bcat-2; Eca40
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009737
ORF Size:	1182 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206133).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009737.1
RefSeq Size:	1780 bp
RefSeq ORF:	1182 bp
Locus ID:	12036
UniProt ID:	O35855
Cytogenetics:	7 29.36 cM
Gene Summary:	Catalyzes the first reaction in the catabolism of the essential branched chain amino acids leucine, isoleucine, and valine. May also function as a transporter of branched chain alpha-keto acids.[UniProtKB/Swiss-Prot Function]



[View online »](#)