

Product datasheet for **MR206159A1V**

Mouse Strada (NM_028126) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Strada (NM_028126) AAV Particle
Tag:	Myc-DDK
Symbol:	Strada
Synonyms:	2610019A05Rik; 6030402H20Rik; A1480680; E130112C09Rik
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR206159 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCTTTTTCTTGCCAATGAGGCGAGCTCAGAGTCCATAGCATCTTTCTCAAACCAGAGATGATGAGCA
GCTTTCTTCCAGAGGGCGGGTGTATGAGCTGCTCACCATCATAGGCAAAGGTTTGAGGACTTGATGAC
AGTGAATTTAGCCAGGTACAAACCAACGGGAGAGTATGTGACAGTACGAAGGATCAACCTAGAAGCTTGT
TCCAATGAGATGGTGACATTCTGCAGGGGAGCTTCATGTCTCTAACTCTTCAGCCATCCCAATATAG
TACCGTATCGAGCCACCTTCATCGCAGATAATGAGCTGTGGGTCGTACATCATTTATGGCGTACGGCTC
TGCAAAGGATCTCATTGGCACACACTTCATGGATGGCATGAATGAAGTGGCGATTGCATACATCCTGCAG
GGGGTGCTAAAGGCCTTAGACTACATCCATCACATGGGCTATGTGCACAGGAGTGTCAAGGCCAGCCACA
TTCTGATCTCCACGGATGGGAAGGTCTACCTGTCTGGTTACGCAGCAACCTTAGCATGATTAGTCACGG
GCAGCGGCAGCGTGAGTCCACGATTTCCCAAATATAGTATCAAGGTCCTGCCTTGCTCAGCCAGAA
GTCCTCCAGCAGAATCTTCAGGGTTACGATGCCAAGTCCGATATCTACAGTGTGGGAATCACAGCATGTG
AACTAGCCAATGGCCATGTCCCCTTCAAAGATATGCCTGCCACTCAGATGCTGCTAGAAAACTTAATGG
CACAGTGCCCTGCCTGCTGGACACCAGCACCATCCCTGCCGAGGAGCTGACCATGAGTCCCTCCGCTCC
ATCGCCAACCCTGGCCTGAATGACAGCCTGGCTGCTGGTAGCCTCCGGCCTTCCAATGGCGACTCGCCTT
CCCACCCCTACCACCGAACCTTCTCGCCTCACTTCCACAACCTTTGTGGAACAGTGCCTTCAGCGCAACCC
TGACGCAAGGCCAAACGCCAGCACCCTGCTGAACCACTCCTTCTCAAGCAGATCAAGCGCGCTGCTTCC
GAGGCTTTGCCTGAGTTGCTCCGCCAGTCACTCCCATCACCCTTCGAGGGCAGCCAGTCTCAGGATC
ACAGTGAATCTTCGGCTGGTGACAAACCTAGAAGACCTGGAGGTGGATGACTGGGAGTTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR206159 protein sequence Red=Cloning site Green=Tags(s) MSFLANEASSESIASFSPKPEMMSSFLPEGGCYELLTIIGKGFEDLMTVNLARYKPTGEYVTVRRINLEAC SNEMVTFLQGELHVSCLFVSHPNIVPYRATFIADNELWVTSFMAYGSAKDLIGTHFMDGMNELAIAYILQ GVLKALDYIHMGYVHRSVKASHILISTDGKVYLSGLRSNLSMISHGQRQRAVHDFPKYSIKVLPWLSPE VLQQNLQGYDAKSDIYSVGITACELANGHVFPKDMPTQMLLEKLNQTVPCLLDTSTIPAEELTMSPSR IANPGLNDSLAAAGSLRPSNGDSPSHPYHRTFSPHFHNFVEQCLQRNPDPARNASTLLNHSFFKQIKRRAS EALPELLRPVTPITNFEQSQSDHSGIFGLVTNLEDLEVDWDF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_028126
ORF Size:	1182 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_028126.1
RefSeq Size:	2432 bp
RefSeq ORF:	1185 bp
Locus ID:	72149
UniProt ID:	Q3UUJ4
Cytogenetics:	11 E1
MW:	43.8 kDa