

Product datasheet for **MR206163L3V**

Esam (NM_027102) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Esam (NM_027102) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Esam
Synonyms:	2310008D05Rik; Esam1; W117m
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_027102
ORF Size:	1182 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206163).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027102.3 , NP_081378.1
RefSeq Size:	1859 bp
RefSeq ORF:	1185 bp
Locus ID:	69524
UniProt ID:	Q925F2
Cytogenetics:	9 A4
Gene Summary:	Can mediate aggregation most likely through a homophilic molecular interaction. [UniProtKB/Swiss-Prot Function]



[View online »](#)