

## Product datasheet for **MR206174L3V**

### Sox11 (NM\_009234) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sox11 (NM_009234) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sox11
Synonyms:	1110038H03Rik; 6230403H02Rik; AI836553; end1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009234
ORF Size:	1185 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206174).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009234.5</a>
RefSeq Size:	7497 bp
RefSeq ORF:	1188 bp
Locus ID:	20666
UniProt ID:	<a href="#">Q7M6Y2</a>
Cytogenetics:	12 A2
Gene Summary:	Transcriptional factor involved in the embryonic neurogenesis. May also have a role in tissue modeling during development.[UniProtKB/Swiss-Prot Function]



[View online »](#)