

Product datasheet for **MR206191**

Sh3glb2 (NM_139302) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sh3glb2 (NM_139302) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Sh3glb2
Synonyms:	SPAS-1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR206191 representing NM_139302 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGACTTCAACATGAAGAAGCTGGCGTCGGACGCGGGCATCTTCTTCACTCGGGCGGTGCAGTTCACAG
AGGAGAAGTTTGGCCAGGCTGAGAAGACGGAGCTTGACGCCCACTTTGAAACCTCTAGCTCGGGCAGA
CAGCACCAAGAACTGGACAGAGCGGATCCTGAGGCAGACCGAGGTGCTGCTGCAGCCCAACCCAGTGCT
CGAGTGGAGGAGTTCTATATGAGAAGCTGGACAGAAAGGTGCCCTCGAGAGTCACCAATGGGGAGCTGC
TGGCTCAGTACATGGCGGAGGCAGCCAGCGAGCTGGGCCCCAGCACTCCCTACGGGAAGACGCTGATCAA
GGTGTCAGAAGCTGAGAAGCGCCTCGGAGCAGCAGAGCGGGATTCATTCACTGCCTCCCTCAGCTTC
CTCACACCCCTGCGGAATTCCTAGAGGGGAGCTGGAAAACGATTTGAAGGAGAGGCGGCTCCTGCAGA
ACCGGCGTCTTGACCTGGATGCCTGCAAAGCCCGGCTAAAGAAGGCCAAGGCAGCCGAAGCCAAAGCCAC
GACGGTGCCTGACTTTCAGGAGACTAGACCTCGTAATTACATTCTATCGGCCAGCGCCTCCGCGCTTTGG
AACGATGAAGTCGACAAGGCTGAGCAGGAGCTTCGAGTGGCGCAGACAGAGTTTGACCGGCAGGCAGAAG
TGACCCGTCTCCTGCTGGAGGGGATCAGCAGTACTCATGTGAACACCTTCGCTGCCTCCATGAGTTCGT
CAAGTCTCAGACAACCTACTATGCACAGTGTACCGCCACATGCTGGATCTGCAGAAACAGCTAGGCAGA
TTTCCAGGCACCTTCGTGGGCACACAGAGCCTGCCTCCCCACCCCTGAGCAGCACCTCACCTACCACCA
CTGCGGCCACCATGCCTGTGGTACCTACTGGGGCTGTCTTGGCCCTCCAGAAGAGGCAGCCCTCTGCCT
GGAGGAGGTGGTCCCCAGCCAGTGGGACTCGAAAGGCCCGGTGCTCTACGACTACGAGGCAGCTGAC
AGCAGTGAGCTGGCCCTGCTGGCTGATGAGCTCATCTGTCTACAGCCTGCCAGGCATGGATCCCGACT
GGCTATTGGTGAGAGAGGCAACAAGAAGGGCAAGGTTCTGTACCTACCTGGAACCTCTCAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR206191 representing NM_139302
Red=Cloning site Green=Tags(s)


MDFNMKKLASDAGIFFTRAVQFTEEFKGQAEKTELDAHFNLLARADSTKNWTERILRQTEVLLQPNPSA
RVEEFLYEKLDRKVP SRVTNGELLAQYMAEAASELGPSTPYGKTLIKVSEAEKRLGAAERDFIHTASLSF
LTPLRNFLEGDWKTISKERRLLQNRRLDL DACKARLKKAKAAEAKATTVPDFQETRPRNYILSASASALW
NDEVDKAEQELRVAQTEFDRQAEVTRLLLEGISSTHVNHLRCLHEFVKSQTTYAQCYRHMLDLQQLGR
FPGTFVGTTEPASPLSSTPTTTAATMPVPTGAVLAPPEEAALCLEEVAPPASGTRKARVLYDYEAAD
SSELALLADELITVYSLPGMDPDWLIGERGNKKGKVPVLYLELLS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_139302

ORF Size:	1185 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_139302.2 , NP_647463.1
RefSeq Size:	1740 bp
RefSeq ORF:	1188 bp
Locus ID:	227700
UniProt ID:	Q8R3V5
Cytogenetics:	2 B
MW:	44.6 kDa