

## Product datasheet for **MR206194L4V**

### Gjc1 (NM\_001159382) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Gjc1 (NM_001159382) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gjc1
Synonyms:	C130009G16Rik; Cnx45; Cx45; Gja-7; Gja7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001159382
ORF Size:	1188 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206194).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001159382.1</a> , <a href="#">NP_001152854.1</a>
RefSeq Size:	1984 bp
RefSeq ORF:	1191 bp
Locus ID:	14615
UniProt ID:	<a href="#">P28229</a>
Cytogenetics:	11 E1
Gene Summary:	One gap junction consists of a cluster of closely packed pairs of transmembrane channels, the connexons, through which materials of low MW diffuse from one cell to a neighboring cell.[UniProtKB/Swiss-Prot Function]



[View online »](#)