

Product datasheet for MR206209L4

H28 (BC024930) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	H28 (BC024930) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	H28
Synonyms:	H-28; H28; H28-1; NS1178
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206209).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	BC024930
ORF Size:	1188 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC024930 , AAH24930
RefSeq Size:	2228 bp
RefSeq ORF:	1190 bp
Locus ID:	15061
Cytogenetics:	3 76.94 cM
Gene Summary:	This gene encodes a minor histocompatibility antigen that is involved in tissue graft rejection and is polymorphic in different mouse strains. The 5' region of this gene, including the translational start codon and the first 293 amino acids, is absent from the C57BL/6J genome. The encoded antigen is not thought to be expressed in C57BL/6J. [provided by RefSeq, Mar 2013]

Product images:



Circular map for MR206209L4