

Protein Sequence:	>MR206216 representing NM_008289 Red=Cloning site Green=Tags(s) MERWPWPSGGAWLLVAARALLQLLRSDRLGRPLLAALALLAALDWLCQRLPPPAALVVLGAGWIALS RLARPPRLPVATRAVLITGCDTGFGKETAKKLDAMGFTVLATVLDLNSPGALELRDLCSPLKLLQMDLT KAEDISRVLEITKAHTASTGLWGLVNNAGLNIVVADVELSPVATFRKCMEVNFFGALELTGKLLPLLRHS RGRIVTVGSPAGDMPYPCLAAYGTSKAAIALLMDFGCELLPWGIKVSIIKPGCFKTDVNTNVLWEKRR QLLLANIPRELLQAYGEDYIEHVHGQFLNSLRMALPDLSPVVDIIDALLAAQPRSRYPGRGLGLMYFI HHYLPEGLRRCFLQNFFINHLLPRALRPGQHGPA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_008289
ORF Size:	1158 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_008289.2 , NP_032315.2
RefSeq Size:	1874 bp
RefSeq ORF:	1161 bp
Locus ID:	15484
UniProt ID:	P51661
Cytogenetics:	8 53.04 cM
MW:	42.6 kDa