

Product datasheet for MR206230L3

Gfra3 (NM_010280) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gfra3 (NM_010280) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Gfra3
Synonyms:	GFRa; GFRalpha3; Y15110
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206230).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_010280
ORF Size:	1191 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_010280.3 , NP_034410.3
RefSeq Size:	1992 bp
RefSeq ORF:	1194 bp
Locus ID:	14587
UniProt ID:	O35118
Cytogenetics:	18 B1
Gene Summary:	This gene encodes a cell surface glycoprotein and member of the glial cell line-derived neurotrophic receptor (GDNFR) family of proteins. The encoded preproprotein is proteolytically processed to generate the mature protein. This protein mediates binding of the ligand artemin to the ret receptor tyrosine kinase, and this interaction may regulate thermal pain and axon regeneration. Homozygous knockout mice for this gene exhibit impaired proliferation of cultured neuroblasts and impaired development of the superior cervical ganglion. [provided by RefSeq, Aug 2015]

Product images:



Circular map for MR206230L3