

Product datasheet for MR206233L3

Vps37a (BC039290) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Vps37a (BC039290) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Vps37a
Synonyms:	2210018P21Rik; 4930592A21Rik; AW261445; D8Ertd531e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206233).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	BC039290
ORF Size:	1191 bp



[View online »](#)

