

## Product datasheet for **MR206241L4V**

### Cyth1 (NM\_001112699) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cyth1 (NM_001112699) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cyth1
Synonyms:	CLM1; CTH-1; CYTIP; Pscd1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001112699
ORF Size:	1194 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206241).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001112699.1</a>
RefSeq Size:	3176 bp
RefSeq ORF:	1194 bp
Locus ID:	19157
UniProt ID:	<a href="#">Q9QX11</a>
Cytogenetics:	11 E2



[View online »](#)

**Gene Summary:**

Promotes guanine-nucleotide exchange on ARF1, ARF5 and ARF6 (PubMed:18042453, PubMed:20080746). Promotes the activation of ARF factors through replacement of GDP with GTP (PubMed:18042453). Plays an important role in membrane trafficking, during junctional remodeling and epithelial polarization, through regulation of ARF6 activity (PubMed:20080746, PubMed:29420262).[UniProtKB/Swiss-Prot Function]