

Product datasheet for **MR206242A1V**

Mouse Chil3 (NM_009892) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Chil3 (NM_009892) AAV Particle
Tag:	Myc-DDK
Symbol:	Chil3
Synonyms:	AI505981; Chi3; Chi3l3; ECF-L; Ym1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR206242 representing NM_009892 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCAAGCTCATTCTTGTCACAGGTCTGGCAATTCTTCTGAACGTACAGCTGGGATCTTCTACCAGC
TGATGTGCTACTATACCAAGTTGGGCTAAGGACAGGCCAATAGAAGGGAGTTTCAAACCTGGTAATATTGA
CCCCTGCCTGTGTAACCTGATCTATGCCTTTGCTGGAATGCAGAATAATGAGATCACTTACACACAT
GAGCAAGACTTGCCTGACTATGAAGCATTGAATGGTCTGAAAGACAAGAACACTGAGCTAAAACTCTCC
TGGCCATTGGAGGATGGAAGTTTGGACCTGCCCGTTTCAGTGCCATGGTCTCTACTCCTCAGAACCGTCA
GATATTCACTTCAGTCAGTTATCAGATTCTTCGTCGAATATAACTTTGATGGCCTCAACCTGGACTGGCAG
TACCCTGGGTCTCGAGGAAGCCCTCCTAAGGACAAACATCTCTTCAGTGTTCTGGTGAAGGAAATGCGTA
AAGCTTTTGGGAAGAATCTGTGGAGAAAGACATTCCAAGGCTGCTACTCACTTCCACAGGAGCAGGAAT
CATTGACGTAATCAAGTCTGGGTACAAGATCCCTGAAGTGTCTCAGTCTCTTGACTATATTCAGGTCATG
ACATATGATCTCCATGATCCTAAGGATGGCTACACTGGAGAAAATAGTCCCTCTATAAATCTCCATATG
ACATTGGAAGAGTGCTGATCTCAATGTGGATTCAATCATTTCTACTGGAAGGACCATGGAGCAGCTTC
TGAGAAGCTCATTGTGGGATTTCCAGCATATGGGCATACCTTTATCCTGAGTGACCTTCTAAGACTGGA
ATTGGTGCCCTACAATTAGTACTGGCCACCAGGAAAGTACACAGATGAATCAGGACTCCTGGCTTACT
ATGAGGTTTGTACATTTCTGAATGAAGGAGCCACTGAGGTCTGGGATGCCCGCCAGGAAGTACCTATGC
CTATCAGGGAATGAGTGGGTTGGTTATGACAATGTCAGGAGCTTCAAGTTGAAGGCTCAGTGGCTCAAG
GACAACAATTTAGGAGGTGCCGTGGTCTGGCCCTGGACATGGATGACTTCAGTGGTTCTTTCTGTCAAC
AGAGACATTTCCCTCTGACATCTACTTTAAAGGAGATCTCAATATACACAGTGAAGTTGAAGGGCCC
TTAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR206242 representing NM_009892 Red=Cloning site Green=Tags(s) MAKLILVTGLAILLNVLGSSYLQMCYYTSWAKDRPIEGSFKPGNIDPCLCTHLYAFAGMQNNEITYTH EQDLRDYEALNGLKDKNTEKTLAIGGWKFGPAPFSAMVSTPQNRQIFIQSVIRFLRQYNFDGLNDWQ YPGSRGSPPKDKHLFSLVKEMRKAEEEEVEKDIPRLLLTSTGAGIIDVIKSGYKIPELSQSLDYIQVM TYDLHDPKDGTYGENSPLYKSPYDIGKSADLNVDISIISYWKDHGAASEKLIVGFPAYGHTFILSDPSKTG IGAPTISTGPPGKYTDESGLLAYVEVCTFLNEGATEVWDAPQEVYAYQGENWGYDNVRSFKLKAQWLK DNNLGGAVVWPLDMDDFSGSFCHQRHFPLTSTLKGDLNIHSASCKGPY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_009892
ORF Size:	1194 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_009892.2
RefSeq Size:	1567 bp
RefSeq ORF:	1197 bp
Locus ID:	12655
UniProt ID:	O35744
Cytogenetics:	3 46.49 cM
MW:	44.9 kDa