

Product datasheet for **MR206245L3V**

Aadac (NM_023383) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Aadac (NM_023383) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Aadac
Synonyms:	5033417E09Rik; Aada; AI265437
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_023383
ORF Size:	1197 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206245).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023383.1 , NP_075872.1
RefSeq Size:	1294 bp
RefSeq ORF:	1197 bp
Locus ID:	67758
UniProt ID:	Q99PG0
Cytogenetics:	3 D



[View online »](#)

Gene Summary:

Displays cellular triglyceride lipase activity in liver, increases the levels of intracellular fatty acids derived from the hydrolysis of newly formed triglyceride stores and plays a role in very low-density lipoprotein assembly (By similarity). Displays serine esterase activity in liver. Deacetylates a variety of arylacetamide substrates, including xenobiotic compounds and procarcinogens, converting them to the primary arylamide compounds and increasing their toxicity.[UniProtKB/Swiss-Prot Function]