

## Product datasheet for **MR206267L3V**

### **Hoxa10 (NM\_008263) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Hoxa10 (NM_008263) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hoxa10
Synonyms:	Hox-1.8; Hoxa-10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008263
ORF Size:	1248 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206267).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008263.3</a> , <a href="#">NP_032289.2</a>
RefSeq Size:	2581 bp
RefSeq ORF:	1251 bp
Locus ID:	15395
UniProt ID:	<a href="#">P31310</a>
Cytogenetics:	6 25.4 cM



[View online »](#)

**Gene Summary:**

In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of a cluster on chromosome 6 and encodes a DNA-binding transcription factor that may regulate gene expression, morphogenesis, and differentiation. More specifically, it may function in fertility, embryo viability, and regulation of hematopoietic lineage commitment. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]