

Product datasheet for **MR206269L3V**

Cyth3 (NM_011182) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cyth3 (NM_011182) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cyth3
Synonyms:	AI648983; ARNO3; CLM3; Grp1; Pscd3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011182
ORF Size:	1197 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206269).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011182.4 , NP_035312.3
RefSeq Size:	3749 bp
RefSeq ORF:	1200 bp
Locus ID:	19159
UniProt ID:	O08967
Cytogenetics:	5 82.36 cM
Gene Summary:	Promotes guanine-nucleotide exchange on ARF1. Promotes the activation of ARF factors through replacement of GDP with GTP (PubMed:18042453). Play a role in the epithelial polarization (PubMed:20080746).[UniProtKB/Swiss-Prot Function]



[View online »](#)