

Product datasheet for **MR206272L3V**

Nolc1 (BC003244) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nolc1 (BC003244) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nolc1
Synonyms:	P130, NOPP130, NOPP140, mKIAA0035
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC003244
ORF Size:	1197 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206272).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC003244
RefSeq Size:	1872 bp
RefSeq ORF:	1199 bp
Locus ID:	70769
Cytogenetics:	19 C3



[View online »](#)

Gene Summary:

Nucleolar protein that acts as a regulator of RNA polymerase I by connecting RNA polymerase I with enzymes responsible for ribosomal processing and modification (By similarity). Required for neural crest specification: following monoubiquitination by the BCR(KBTBD8) complex, associates with TCOF1 and acts as a platform to connect RNA polymerase I with enzymes responsible for ribosomal processing and modification, leading to remodel the translational program of differentiating cells in favor of neural crest specification (By similarity). Involved in nucleologenesis, possibly by playing a role in the maintenance of the fundamental structure of the fibrillar center and dense fibrillar component in the nucleolus (PubMed:11424213). It has intrinsic GTPase and ATPase activities (By similarity). [UniProtKB/Swiss-Prot Function]