

Product datasheet for **MR206281L4V**

Klhl7 (NM_001161800) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Klhl7 (NM_001161800) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Klhl7
Synonyms:	2700038B03Rik; D5Ertd363e; SBBI26
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001161800
ORF Size:	1203 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206281).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001161800.1
RefSeq Size:	3329 bp
RefSeq ORF:	1422 bp
Locus ID:	52323
Cytogenetics:	5 10.67 cM
Gene Summary:	Substrate-specific adapter of a BCR (BTB-CUL3-RBX1) E3 ubiquitin ligase complex. The BCR(KLHL7) complex acts by mediating ubiquitination and subsequent degradation of substrate proteins. Probably mediates 'Lys-48'-linked ubiquitination (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)