

Product datasheet for **MR206290L4V**

Fbxl14 (NM_133940) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Symbol: Fbxl14

Synonyms: AW322056; Fbx14l

Mammalian Cell Selection: Puromycin

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_133940

ORF Size: 1200 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR206290).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_133940.1](#)

RefSeq Size: 4219 bp

RefSeq ORF: 1203 bp

Locus ID: 101358

UniProt ID: [Q8BID8](#)

Cytogenetics: 6 F1



Gene Summary:

Substrate-recognition component of some (SKP1-CUL1-F-box protein)-type E3 ubiquitin-protein ligase complexes. The SCF(FBXL14) complex acts by mediating ubiquitination and subsequent degradation of SNAIL (By similarity).[UniProtKB/Swiss-Prot Function]