

Product datasheet for **MR206299L3V**

Kcnab1 (NM_010597) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kcnab1 (NM_010597) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kcnab1
Synonyms:	Akr8a8; Kvbeta1.1; mKv(beta)1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010597
ORF Size:	1206 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206299).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010597.2
RefSeq Size:	3267 bp
RefSeq ORF:	1206 bp
Locus ID:	16497
UniProt ID:	P63143
Cytogenetics:	3 30.15 cM



[View online »](#)

Gene Summary:

Cytoplasmic potassium channel subunit that modulates the characteristics of the channel-forming alpha-subunits (PubMed:10454353). Modulates action potentials via its effect on the pore-forming alpha subunits (PubMed:10454353). Promotes expression of the pore-forming alpha subunits at the cell membrane, and thereby increases channel activity (PubMed:8824288). Mediates closure of delayed rectifier potassium channels by physically obstructing the pore via its N-terminal domain and increases the speed of channel closure for other family members (By similarity). Promotes the closure of KCNA1, KCNA2 and KCNA5 channels (By similarity). Accelerates KCNA4 channel closure (By similarity). Accelerates the closure of heteromeric channels formed by KCNA1 and KCNA4 (By similarity). Accelerates the closure of heteromeric channels formed by KCNA2, KCNA5 and KCNA6 (By similarity). Enhances KCNB1 and KCNB2 channel activity (PubMed:8824288). Binds NADPH; this is required for efficient down-regulation of potassium channel activity (By similarity). Has NADPH-dependent aldo-ketoreductase activity (By similarity). Oxidation of the bound NADPH strongly decreases N-type inactivation of potassium channel activity (By similarity). [UniProtKB/Swiss-Prot Function]