

Product datasheet for **MR206300L4V**

C730036D15Rik (BC010829) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C730036D15Rik (BC010829) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	C730036D15Rik
Synonyms:	C730036D15Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC010829
ORF Size:	1203 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206300).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC010829 , AAH10829
RefSeq Size:	1442 bp
RefSeq ORF:	1205 bp
Locus ID:	209186
Cytogenetics:	4 B1
Gene Summary:	Acytransferase which efficiently conjugates very long-chain and long-chain fatty acids to taurine. Shows no conjugation activity in the presence of glycine (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)