

Product datasheet for **MR206334L4V**

Lipk (BC055815) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lipk (BC055815) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lipk
Synonyms:	9930115F20; BC055815; LipI2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC055815
ORF Size:	1209 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206334).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC055815 , AAH55815
RefSeq Size:	2136 bp
RefSeq ORF:	1211 bp
Locus ID:	240633
Cytogenetics:	19 C1
Gene Summary:	Plays a highly specific role in the last step of keratinocyte differentiation. May have an essential function in lipid metabolism of the most differentiated epidermal layers (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)