

Product datasheet for **MR206365L4V**

Poc1a (NM_027354) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Poc1a (NM_027354) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Poc1a
Synonyms:	2510040D07Rik; cha; Wdr51a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_027354
ORF Size:	1218 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206365).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027354.1 , NP_081630.1
RefSeq Size:	1707 bp
RefSeq ORF:	1218 bp
Locus ID:	70235
UniProt ID:	Q8JZX3
Cytogenetics:	9 F1


[View online »](#)

Gene Summary:

Plays an important role in centriole assembly and/or stability and ciliogenesis. Involved in early steps of centriole duplication, as well as in the later steps of centriole length control. Acts in concert with POC1B to ensure centriole integrity and proper mitotic spindle formation (By similarity).[UniProtKB/Swiss-Prot Function]