

Product datasheet for **MR206367L4V**

Rbm17 (NM_152824) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rbm17 (NM_152824) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rbm17
Synonyms:	2700027J02Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_152824
ORF Size:	1218 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206367).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_152824.1 , NP_690037.1
RefSeq Size:	1583 bp
RefSeq ORF:	1218 bp
Locus ID:	76938
UniProt ID:	Q8JZX4
Cytogenetics:	2 A1
Gene Summary:	Splice factor that binds to the single-stranded 3'AG at the exon/intron border and promotes its utilization in the second catalytic step. Involved in the regulation of alternative splicing and the utilization of cryptic splice sites (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)