

Product datasheet for **MR206406L4V**

Snapc3 (NM_029949) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Snapc3 (NM_029949) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Snapc3
Synonyms:	1810020H02Rik; 4930558A07Rik; 5031401C21Rik; AI414457; E030018J20Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_029949
ORF Size:	1224 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206406).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029949.1
RefSeq Size:	1413 bp
RefSeq ORF:	1224 bp
Locus ID:	77634
UniProt ID:	Q9D2C9
Cytogenetics:	4 C3


[View online »](#)

Gene Summary:

Part of the SNAPc complex required for the transcription of both RNA polymerase II and III small-nuclear RNA genes. Binds to the proximal sequence element (PSE), a non-TATA-box basal promoter element common to these 2 types of genes. Recruits TBP and BRF2 to the U6 snRNA TATA box (By similarity).[UniProtKB/Swiss-Prot Function]