

## Product datasheet for **MR206423L3V**

### Cas1 (BC035954) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cas1 (BC035954) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cas1
Synonyms:	2410019P08Rik; AV096627; BC040081; Cst; D4Ertd432e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC035954
ORF Size:	1227 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206423).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC035954</a> , <a href="#">AAH35954</a>
RefSeq Size:	2091 bp
RefSeq ORF:	1229 bp
Locus ID:	69743
Cytogenetics:	4 78.87 cM
Gene Summary:	Transcription factor involved in vascular assembly and morphogenesis through direct transcriptional regulation of EGFL7.[UniProtKB/Swiss-Prot Function]



[View online »](#)