

Product datasheet for **MR206430L3V**

Elk3 (NM_013508) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Elk3 (NM_013508) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Elk3
Synonyms:	D430049E23Rik; Erp; Etrp; Net; Sap-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013508
ORF Size:	1230 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206430).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013508.1 , NP_038536.1
RefSeq Size:	4234 bp
RefSeq ORF:	1230 bp
Locus ID:	13713
UniProt ID:	P41971
Cytogenetics:	10 48.04 cM



[View online »](#)

Gene Summary:

May be a negative regulator of transcription, but can activate transcription when coexpressed with Ras, Src or Mos. Forms a ternary complex with the serum response factor and the ETS and SRF motifs of the Fos serum response element.[UniProtKB/Swiss-Prot Function]