

Product datasheet for **MR206440L4V**

Fam111a (NM_026640) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fam111a (NM_026640) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fam111a
Synonyms:	4632417K18Rik; AW413625
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026640
ORF Size:	1233 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206440).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026640.2 , NP_080916.1
RefSeq Size:	3199 bp
RefSeq ORF:	1842 bp
Locus ID:	107373
UniProt ID:	Q9D2L9
Cytogenetics:	19
Gene Summary:	Chromatin-associated protein required for PCNA loading on replication sites. Promotes S-phase entry and DNA synthesis (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)