

## Product datasheet for **MR206443L4V**

### Dpep1 (NM\_007876) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Dpep1 (NM_007876) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dpep1
Synonyms:	AI327012; MBD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007876
ORF Size:	1230 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206443).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_007876.1</a> , <a href="#">NP_031902.1</a>
RefSeq Size:	2249 bp
RefSeq ORF:	1233 bp
Locus ID:	13479
UniProt ID:	<a href="#">P31428</a>
Cytogenetics:	8 72.08 cM
Gene Summary:	Hydrolyzes a wide range of dipeptides. Implicated in the renal metabolism of glutathione and its conjugates. Converts leukotriene D4 to leukotriene E4; it may play an important role in the regulation of leukotriene activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)