

Product datasheet for **MR206474L4V**

Dnaja2 (NM_019794) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dnaja2 (NM_019794) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dnaja2
Synonyms:	1500017M13Rik; 2010206B19Rik; DNAJ; Dnaj3; DNJ3; HIRIP4; mDj3; PRO3015
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019794
ORF Size:	1239 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206474).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019794.1 , NP_062768.1
RefSeq Size:	2839 bp
RefSeq ORF:	1239 bp
Locus ID:	56445
UniProt ID:	Q9QYJ0
Cytogenetics:	8 C3
Gene Summary:	Co-chaperone of Hsc70. Stimulates ATP hydrolysis and the folding of unfolded proteins mediated by HSPA1A/B (in vitro).[UniProtKB/Swiss-Prot Function]



[View online »](#)