

Product datasheet for **MR206488**

Suv39h1 (NM_011514) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Suv39h1 (NM_011514) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Suv39h1
Synonyms:	AI852103; AL022883; DXHXS7466e; H3-K9-HMTase 1; KMT1A; miS6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR206488 representing NM_011514 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGAAAATTTAAAAGGTTGCAGTGTGTGCTGTAATCTTCTTGAATCAACTGCAGGACCTGTGCC
GACTAGCCAAGCTTTCTGTCTGCCCTGGTGTCTAAGAAGAATCTGTATGACTTTGAAGTTGAATA
CCTGTGTGATTATAAGAAGATCCGTGAGCAGGAGTATTACCTGGTTAAGTGGCGTGGGTATCCCGACTCA
GAAAACACCTGGGAGCCACGGCAGAATCTAAAATGTATACGAGTCTTAAGCAGTTCACAAGGACTTAG
AAAGAGAGCTTGTCCGACGACACCGCCGGTCAAAGCCACCCAGGCATCTGGACCCAAACCTAGCCAATTA
CCTGGTGCAGAAGGCCAAGCAGAGGGCGGCACTTCAGCGTTGGGAACAAGAGCTCAATGCCAAGCGCAGC
CACCTGGGGCGGATCACCGTGGAGAATGAGGTAGACCTGGATGGCCCTCAAGGTCTTTGTCTATATCA
ATGAGTATCGAGTTGGTGAAGGCATCACCTCAACCAGGTAGCTGTTGGCTGTGAGTGCCAGGACTGTCT
GTTGGCACCCACTGGAGGCTGTTGCCCTGGAGCATCCCTGCACAAGTTGCCTACAATGACCAAGGCCAG
GTGCGACTGAAAGCTGGGAGCCATCTACGAGTGCACTCCCGCTGTTGCTGTGGCTATGACTGCCCAA
ACCGTGTAGTCCAGAAAGGCATCCGCTACGATCTCTGCATCTCCGCACTAATGATGGCCGAGGCTGGGG
TGTCCGCAGCTGAAAAGATCCGAAAAATAGCTTTGTTATGGAGTATGTGGGAGAGATTATTACCTCA
GAGGAGCAGAGCGGAGGGGCCAGATCTACGACCGCCAGGGCCACCTACCTTTTGACCTGGACTACG
TGGAAAGACGTATATACCGTGGATGCCGCTTATTATGGCAACATCTCTATTTTGTCAACCATAGTTGTGA
TCCCAACCTGCAGGTGTACAACGTATTCATAGACAACCTTGATGAGCGACTACCCCGCATCGCATTCTTT
GCCACAAGAACCATCTGGGCGGGCGAGGAGCTCACCTTTGATTACAACATGCAAGTGGACCCCGTGGACA
TGGAGAGTACCCGAATGGACTCCAACCTTTGGCCTGGCTGGGCTCCCCGGCTCCCCAAGAAACGAGTCCG
TATTGAATGCAAATGTGGGACAACGGCTTGCCGAAAATACCTCTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR206488 representing NM_011514
 Red=Cloning site Green=Tags(s)

MAENLKGCSVCKSSWNQLQDLRLAKLSCPALGVSKKNLYDFEVEYLCDYKKIREQEYYLVKWRGYPDS
 ENTWEPRQNLKCIRVLKQFHKLDERELVRRHRRSKPPRHLDPNLANYLQKAKQRRALQRWEQELNAKRS
 HLGRLITVENEVDLDGPPRSFVYINEYRVGEGITLNQVAVGCECQDCLLAPTGGCCPGASLHKFAYNDAQGQ
 VRLKAGQPIYECNSRCCCGYDCPNRVVQKGIKYDLICIFRTNDGRGWGVRTLEKIRKNSFVMEYVGEIITS
 EEAERRGQIYDRQGATYLFDLDYVEDVYTVDAAYYGNISHFVNHSCDPNLQVYVNFIDNLDLRLPRIAFF
 ATRTIWAGEELTFDYNMQVDPVDMESTRMSNFGLAGLPGSPKKRVRIECKGTTACRKLYL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_011514

ORF Size: 1236 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_011514.2](#)

RefSeq Size: 2820 bp

RefSeq ORF: 1239 bp

Locus ID: 20937


UniProt ID: [O54864](#)

Cytogenetics: X 3.64 cM

MW: 48.2 kDa

Gene Summary: Histone methyltransferase that specifically trimethylates 'Lys-9' of histone H3 using monomethylated H3 'Lys-9' as substrate. H3 'Lys-9' trimethylation represents a specific tag for epigenetic transcriptional repression by recruiting HP1 (CBX1, CBX3 and/or CBX5) proteins to methylated histones. Mainly functions in heterochromatin regions, thereby playing a central role in the establishment of constitutive heterochromatin at pericentric and telomere regions. H3 'Lys-9' trimethylation is also required to direct DNA methylation at pericentric repeats. SUV39H1 is targeted to histone H3 via its interaction with RB1 and is involved in many processes, such as repression of MYOD1-stimulated differentiation, regulation of the control switch for exiting the cell cycle and entering differentiation, repression by the PML-RARA fusion protein, BMP-induced repression, repression of switch recombination to IgA and regulation of telomere length. Component of the eNoSC (energy-dependent nucleolar silencing) complex, a complex that mediates silencing of rDNA in response to intracellular energy status and acts by recruiting histone-modifying enzymes. The eNoSC complex is able to sense the energy status of cell: upon glucose starvation, elevation of NAD(+)/NADP(+) ratio activates SIRT1, leading to histone H3 deacetylation followed by dimethylation of H3 at 'Lys-9' (H3K9me2) by SUV39H1 and the formation of silent chromatin in the rDNA locus. Recruited by the PER complex to the E-box elements of the circadian target genes such as PER2 itself or PER1, contributes to the conversion of local chromatin to a heterochromatin-like repressive state through H3 'Lys-9' trimethylation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR206488