

Product datasheet for MR206535

Kcnk2 (BC062094) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Kcnk2 (BC062094) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Kcnk2
Synonyms:	A430027H14Rik; AI848635; TREK-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR206535 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGTGCAGTGGCGGCCCTGACTTGCTGGATCCCAAGTCTGCTGCTCAGAACTCCAAACCGAGGCTCT
CATTCTCCTCAAACCCACCGTCTTGCCTCCCGGTGGAGAGTACTCGGCCATTAATGTTATGAAATG
GAAGACAGTCTCCACGATTTTCTGGTGGTCGCTCTACCTGATCATCGGAGCCACGGTGTCAAGGCA
TTGGAGCAGCCTCAGGAGATTTCCAGAGGACCACCATTGTGATCCAGAAGCAGACCTTCATAGCCAGC
ATGCCTGCGTCACTCCACCGAGCTGGACGAACATCCAGCAAATAGTGGCAGCAATAAACGCAGGGAT
TATCCCCTTAGGAAACAGCTCCAATCAAGTTAGTCACTGGGACCTCGGAAGCTCTTTCTTCTTCTGCTGGT
ACTGTTATCACAACCATAGGATTTGAAACATCTCCCCACGAACCTGAAGGTGGAAAAATATTCTGCATCA
TCTATGCCTTGGTGGAAATCCCCTCTTGGCTTCTACTGGCTGGGGTGGTGTATCAGCTAGGAACTAT
ATTTGGAAAAGGAATTGCCAAAGTGAAGACACATTTATTAAGTGAATGTTAGTCAGACGAAGATTCGT
ATCATCTCCACCATCATCTTCATCCTGTTGGCTGTGCTCTTTGTGGCTCTCCCTGCGGTATATTCA
AGCACATAGAAGGCTGGAGCGCCTGGACGCTATCTATTTGTGGTTATCACTCTGACGACCATTGGATT
TGGAGACTACGTGGCAGGTGGATCAGACATTGAATATCTGGACTTCTACAAGCCTGTGGTGTGGTCTGG
ATCCTCGTTGGGCTGGCCTACTTTGCAGCTGTTCTGAGCATGATTGGGGACTGGCTACGGGTGATCTCA
AGAAGACGAAGGAAGAGGTGGGAGAGTTAGAGCGCATGCCGCTGAGTGGACAGCCAATGTACGGCCGA
GTTCAAGGAAACGAGGAGCGGCTGAGCGTGGAGATCTACGACAAGTCCAGCGTGCCACATCCGTGAAG
CGGAAGCTCTCCGAGAGCTGGCGGGCAACCACAACCAGGAAGTACTCCGTGTAGGAGGACCCTGTCTG
TGAACCACCTGACCAGCGAGAGGAAAGTCTGCCTCCCTTGTGAAGGCTGAGAGCATCTATCTGAACGG
TCTGACACCACACTGTGCTGGTGGAGACATAGCTGTCATTGAGAATGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR206535 protein sequence
Red=Cloning site Green=Tags(s)

MGAVAAPDLLDPKSAQNSKPRLSFSSKPTVLASRVESDSAINVMKWTVSTIFLVVLYLIIGATVFKA
 LEQPQEISQRTTIVIQKQTFIAQHACVNSTELDELIQQIVAAINAGIPLGNSSNQVSHWDLGSSFFFAG
 TVITTIIGFGNISPRTEGGKIFCIIYALLGIPLFGFLLAGVGDQLGTIFGKGIKVEDTFIKWNVSQTKIR
 IISTIIIFILFGCVLFVALPAVIFKHIEGWSALDAIYFVVITLTTIGFGDYVAGGSDIEYLDYFKPVVWF
 ILVGLAYFAAVLSMIGDWLRVISKKTKEEVGEFRAHAAEWTANVTAEFKETRRRLSVEIYDKFQRATSVK
 RKLSAELAGHNQELTPCRRTLSVNHLTSEREVLPPLLKAESIYLNGLTPHCAGEDIAVIENMK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: BC062094

ORF Size: 1242 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC062094](#), [AAH62094](#)

RefSeq Size: 3228 bp

RefSeq ORF: 1244 bp


Locus ID: 16526

Cytogenetics: 1 H6

MW: 45.6 kDa

Gene Summary: Ion channel that contributes to passive transmembrane potassium transport. Reversibly converts between a voltage-insensitive potassium leak channel and a voltage-dependent outward rectifying potassium channel in a phosphorylation-dependent manner. In astrocytes, forms mostly heterodimeric potassium channels with KCNK1, with only a minor proportion of functional channels containing homodimeric KCNK2 (PubMed:24496152). In astrocytes, the heterodimer formed by KCNK1 and KCNK2 is required for rapid glutamate release in response to activation of G-protein coupled receptors, such as F2R and CNR1 (PubMed:24496152).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR206535