

Product datasheet for **MR206537L4V**

Scrn1 (NM_027268) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Scrn1 (NM_027268) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Scrn1
Synonyms:	2810019K23Rik; 6330535A03Rik; AI852905; mKIAA0193; SES1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_027268
ORF Size:	1245 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206537).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027268.1
RefSeq Size:	1777 bp
RefSeq ORF:	1245 bp
Locus ID:	69938
UniProt ID:	Q9CZC8
Cytogenetics:	6 B3
Gene Summary:	Regulates exocytosis in mast cells. Increases both the extent of secretion and the sensitivity of mast cells to stimulation with calcium (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)