

Product datasheet for **MR206577L3V**

BC023179 (BC108353) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	BC023179
Synonyms:	MGC37070
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC108353
ORF Size:	1248 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR206577).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	BC108353 , AA108354
RefSeq Size:	1703 bp
RefSeq ORF:	1250 bp
Locus ID:	232855
Cytogenetics:	7 A1

