

Product datasheet for **MR206585L4V**

Stk26 (NM_133729) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Stk26 (NM_133729) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Stk26
Synonyms:	2610018G03Rik; Mst4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133729
ORF Size:	1251 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206585).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133729.1 , NP_598490.1
RefSeq Size:	2016 bp
RefSeq ORF:	1251 bp
Locus ID:	70415
UniProt ID:	Q99JT2
Cytogenetics:	X A5
Gene Summary:	Mediator of cell growth. Modulates apoptosis. In association with STK24 negatively regulates Golgi reorientation in polarized cell migration upon RHO activation.[UniProtKB/Swiss-Prot Function]



[View online »](#)