

## Product datasheet for **MR206588L3V**

### Hat1 (NM\_026115) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hat1 (NM_026115) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hat1
Synonyms:	2410071B14Rik; AA536933; Hat-1; KAT1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026115
ORF Size:	1251 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206588).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026115.2</a> , <a href="#">NP_080391.2</a>
RefSeq Size:	1576 bp
RefSeq ORF:	1251 bp
Locus ID:	107435
UniProt ID:	<a href="#">Q8BY71</a>
Cytogenetics:	2 C2



[View online »](#)

**Gene Summary:**

Acetylates soluble but not nucleosomal histone H4 at 'Lys-5' (H4K5ac) and 'Lys-12' (H4K12ac) and, to a lesser extent, acetylates histone H2A at 'Lys-5' (H2AK5ac). Has intrinsic substrate specificity that modifies lysine in recognition sequence GXGKXG. May be involved in nucleosome assembly during DNA replication and repair as part of the histone H3.1 and H3.3 complexes.[UniProtKB/Swiss-Prot Function]