

Product datasheet for **MR206605L3V**

Prim1 (NM_008921) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Prim1 (NM_008921) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Prim1
Synonyms:	AI324982
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008921
ORF Size:	1251 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206605).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008921.2 , NP_032947.1
RefSeq Size:	1533 bp
RefSeq ORF:	1254 bp
Locus ID:	19075
UniProt ID:	P20664
Cytogenetics:	10 D3
Gene Summary:	DNA primase is the polymerase that synthesizes small RNA primers for the Okazaki fragments made during discontinuous DNA replication.[UniProtKB/Swiss-Prot Function]



[View online »](#)