

Product datasheet for **MR206632L4V**

Adrb2 (NM_007420) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adrb2 (NM_007420) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Adrb2
Synonyms:	Adrb; Adrb-2; Bad; Badm; Gpc; Gpcr7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007420
ORF Size:	1257 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206632).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007420.2
RefSeq Size:	2269 bp
RefSeq ORF:	1257 bp
Locus ID:	11555
UniProt ID:	P18762
Cytogenetics:	18 35.1 cM



[View online »](#)

Gene Summary:

This intronless gene belongs to the G-protein-coupled receptor superfamily, which includes transmembrane proteins that play a role in signal transduction across biological membranes resulting in a variety of physiological responses. The encoded protein is a beta-2 adrenergic receptor which is activated by catecholamine ligands such as adrenaline and epinephrine. The protein participates in the classical signaling pathway involving G protein, adenylyl cyclase, cAMP (3'-5'-cyclic adenosine monophosphate) and protein kinase A (PKA). In humans, this gene is implicated in susceptibility to nocturnal asthma, obesity and type 2 diabetes. [provided by RefSeq, Apr 2013]