

Product datasheet for **MR206678L3V**

Uxs1 (NM_026430) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Uxs1 (NM_026430) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Uxs1
Synonyms:	1600025113Rik; AI451869; AI649125; AW550562
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026430
ORF Size:	1263 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206678).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026430.1 , NP_080706.1
RefSeq Size:	1263 bp
RefSeq ORF:	1263 bp
Locus ID:	67883
UniProt ID:	Q91XL3
Cytogenetics:	1 C1.1
Gene Summary:	Catalyzes the NAD-dependent decarboxylation of UDP-glucuronic acid to UDP-xylose. Necessary for the biosynthesis of the core tetrasaccharide in glycosaminoglycan biosynthesis (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)