

Product datasheet for **MR206687L3V**

Fam110c (NM_027828) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fam110c (NM_027828) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fam110c
Synonyms:	9030611O19Rik; AU019489
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_027828
ORF Size:	1266 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206687).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027828.2 , NP_082104.2
RefSeq Size:	2684 bp
RefSeq ORF:	1266 bp
Locus ID:	104943
UniProt ID:	Q8VE94
Cytogenetics:	12 A2
Gene Summary:	May play a role in microtubule organization. May play a role in cell spreading and cell migration of epithelial cells; the function may involve the AKT1 signaling pathway. [UniProtKB/Swiss-Prot Function]



[View online »](#)