

Product datasheet for **MR206704L4V**

Acadm (NM_007382) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Acadm (NM_007382) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Acadm
Synonyms:	AU018656; MCA; MCAD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007382
ORF Size:	1266 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206704).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007382.1 , NP_031408.1
RefSeq Size:	2062 bp
RefSeq ORF:	1266 bp
Locus ID:	11364
UniProt ID:	P45952
Cytogenetics:	3 78.77 cM



[View online »](#)

Gene Summary:

This gene encodes a homotetrameric mitochondrial flavoprotein and is a member of the acyl-CoA dehydrogenase family. Members of this family catalyze the first step of fatty acid beta-oxidation, forming a C2-C3 trans-double bond in a FAD-dependent reaction. As beta-oxidation cycles through its four steps, each member of the acyl-CoA dehydrogenase family works at an optimum fatty acid chain-length. This enzyme has its optimum length between C6- and C12-acylCoA. In mice, deficiency of this gene can cause neonatal mortality as well as fasting and cold intolerance. This gene has multiple, intronless pseudogenes. [provided by RefSeq, Nov 2012]