

Product datasheet for **MR206713LV**

9630058J23Rik (BC034107) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	9630058J23Rik (BC034107) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	9630058J23Rik
Synonyms:	MGC31547
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC034107
ORF Size:	1266 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206713).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC034107.1
RefSeq Size:	3579 bp
RefSeq ORF:	1268 bp
Locus ID:	226744
Cytogenetics:	1 H4
Gene Summary:	Required for targeting of connexins to the plasma membrane.[UniProtKB/Swiss-Prot Function]


[View online »](#)