

## Product datasheet for **MR206730L4V**

### Wdr12 (NM\_021312) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Wdr12 (NM_021312) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Wdr12
Synonyms:	4933402C23Rik; Ytm1; Ytm1p
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021312
ORF Size:	1272 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206730).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021312.4</a>
RefSeq Size:	2946 bp
RefSeq ORF:	1272 bp
Locus ID:	57750
UniProt ID:	<a href="#">Q9JJA4</a>
Cytogenetics:	1 C2
Gene Summary:	Component of the PeBoW complex, which is required for maturation of 28S and 5.8S ribosomal RNAs and formation of the 60S ribosome.[UniProtKB/Swiss-Prot Function]



[View online »](#)