

Product datasheet for **MR206766**

Mto1 (BC063256) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mto1 (BC063256) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Mto1
Synonyms:	2310039H01Rik; 5730419A02Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR206766 representing BC063256
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC


ATGTTCTTCTCCGAGGCCGTGGCCATTGGGCCGAGCGTCCCTCGGCAGACGGCTGCCTTTGCGTCGCC
TCCGCAGCGACAGCGCGCGCCCTGCACCCGCACTTCGACGTGGTAGTTATCGGAGGAGGGCACGCGGG
GACCGAGGCGGCCGCGCGCCGCTCGCTGTGGCTCGCGGACCCTGCTTCTCACGCACCGCGTGGACACG
ATCGGTGAGATGTCTGTAATCCTTCTTTGGTGGCATCGGAAAGGGGCATTTGATGAGAGAAGTGGATG
CCTTGGATGGTCTGTGTTCCCGGATCTGTGACCACTCTGGTATACATTACAAAGTGCTAAACCGCGTAA
GGGCCCTGCTGTGTGGGACTGAGAGCACAGATTGACCGAACTTTATAACAGAACATGCAGAAAGAA
ATCCTGAGCACGCCCTACTCACTGTTCAGAAGGGAGCTGTAGAAGATCTAGTTCTTGCAGAACGAGAGC
CTGGGTACCTGGGAAGTCCAGGGTCAGAGGTGTTGTCTTGGCAGATGGAAGTACAATTTATGCAGAGAG
TGTAATTCTGACTACTGGGACATTTCTGAGAGGCATGATCATAATTGGATTGGAGATGCATCCTGCAGGA
CGTTTGGGGGATCAGCCCTCCATAGGGTTGGCTCAGACTCTGGAGAAGTTGGGATTTATGGTAGGAAGAC
TAAAGACTGGAATCCACCCGACTCGCTAAAGAGTCCATTAATTTCTCCATTCTAAACAAGCACACACC
AGATGATCCATCCATACCATTCAGCTTTCTCAGTGACAGCGTGTGGATCAAGCCGGAAGACCAACTGCCA
TGTTACTTGACTCACACCAACCCCGAGTGGATGCGATTGTCTTGAGAACCTTCACCTGAATAGCCACG
TTCAGGAAACTACAAAAGGGCCCGATACTGTCCCTCTATTGAATCAAAGGTTTTACGTTTTCCAAATCG
TCTGCATCAGGTTTGGTTGGAGCCTGAAGGAATGGATTCCGACCTCATCTACCCGCAAGGGTTATCTGTT
ACACTGCCAGCTGAGTTACAGGAGAAAAATGATAACATGCATCAGAGGCTTGGAGAAAGCCAAAATGGTTC
ACCCAGGCTATGGTGTGCAGTATGATTATTTAGACCCCGTCAGATCAGTCCTTCCCTCGAGACACACTT
GGTTCAGCGGCTGTTCTTTGCTGGACAGATAAATGGAACCACTGGCTATGAGGAAGCTGCAGCTCAGTCC
TCTTCTAGGGGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:
Cloning Scheme:

Sgfl-MluI

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN:


BC063256

ORF Size:

1272 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC063256
RefSeq Size:	5221 bp
RefSeq ORF:	1274 bp
Locus ID:	68291
Cytogenetics:	9 E1
MW:	191.4 kDa
Gene Summary:	Involved in the 5-carboxymethylaminomethyl modification (mnm(5)s(2)U34) of the wobble uridine base in mitochondrial tRNAs.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR206766