

## Product datasheet for **MR206774L4V**

### **Tmem184a (NM\_001161548) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Tmem184a (NM_001161548) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmem184a
Synonyms:	Sdmg1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001161548
ORF Size:	1278 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206774).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001161548.1</a> , <a href="#">NP_001155020.1</a>
RefSeq Size:	2039 bp
RefSeq ORF:	1350 bp
Locus ID:	231832
UniProt ID:	<a href="#">Q3UFJ6</a>
Cytogenetics:	5 G2
Gene Summary:	Acts as a heparin receptor in vascular cells (By similarity). May be involved in vesicle transport in exocrine cells and Sertoli cells (PubMed:18321981, PubMed:19097053).[UniProtKB/Swiss-Prot Function]



[View online »](#)