

Product datasheet for **MR206789L4V**

Mfge8 (NM_001045489) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mfge8 (NM_001045489) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mfge8
Synonyms:	AA408458; AI325141; lactadherin; MFG-E8; Mfgm; MP47; P47; SED1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001045489
ORF Size:	1281 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206789).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001045489.1 , NP_001038954.1
RefSeq Size:	2032 bp
RefSeq ORF:	1281 bp
Locus ID:	17304
Cytogenetics:	7 44.9 cM



[View online »](#)

Gene Summary:

Contributes to phagocytic removal of apoptotic cells in many tissues. Specific ligand for the alpha-v/beta-3 and alpha-v/beta-5 receptors. Also binds to phosphatidylserine-enriched cell surfaces in a receptor-independent manner. Zona pellucida-binding protein which may play a role in gamete interaction (By similarity). Plays an important role in the maintenance of intestinal epithelial homeostasis and the promotion of mucosal healing. Promotes VEGF-dependent neovascularization.[UniProtKB/Swiss-Prot Function]