

## Product datasheet for **MR206815L4V**

### Ednra (NM\_010332) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ednra (NM_010332) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ednra
Synonyms:	AEA001; ET-AR; ETa; Gpcr10; Mhdaaea1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010332
ORF Size:	1284 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206815).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_010332.2</a> , <a href="#">NP_034462.1</a>
RefSeq Size:	3643 bp
RefSeq ORF:	1284 bp
Locus ID:	13617
UniProt ID:	<a href="#">Q61614</a>
Cytogenetics:	8 C1
Gene Summary:	Receptor for endothelin-1. Mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system. The rank order of binding affinities for ET-A is: ET1 > ET2 >> ET3 (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)