

Product datasheet for **MR206819L4V**

Igsf11 (NM_170599) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Igsf11 (NM_170599) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Igsf11
Synonyms:	1700025L02Rik; BT-IgSF; Igsf13
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_170599
ORF Size:	1287 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206819).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_170599.2 , NP_733548.2
RefSeq Size:	3430 bp
RefSeq ORF:	1287 bp
Locus ID:	207683
UniProt ID:	P0C673
Cytogenetics:	16 B4
Gene Summary:	Functions as a cell adhesion molecule through homophilic interaction. Stimulates cell growth (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)