

## Product datasheet for **MR206840L4V**

### Elk1 (NM\_007922) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Elk1 (NM_007922) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Elk1
Synonyms:	Elk-; Elk-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007922
ORF Size:	1290 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206840).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_007922.4</a>
RefSeq Size:	3548 bp
RefSeq ORF:	1290 bp
Locus ID:	13712
UniProt ID:	<a href="#">P41969</a>
Cytogenetics:	X 16.45 cM



[View online »](#)

**Gene Summary:**

This gene is a member of the Ets family of transcription factors and of the ternary complex factor (TCF) subfamily. Proteins of the TCF subfamily form a ternary complex by binding to the the serum response factor and the serum response element in the promoter of the c-fos proto-oncogene. The protein encoded by this gene is a nuclear target for the ras-raf-MAPK signaling cascade. This gene may produce multiple isoforms by the use of alternative translational start codons. [provided by RefSeq, Mar 2012]