

Product datasheet for **MR206873L3V**

Cstf1 (NM_024199) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cstf1 (NM_024199) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cstf1
Synonyms:	1700057K18Rik; AI788832
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024199
ORF Size:	1296 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206873).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024199.2
RefSeq Size:	1820 bp
RefSeq ORF:	1296 bp
Locus ID:	67337
UniProt ID:	Q99LC2
Cytogenetics:	2 94.84 cM
Gene Summary:	One of the multiple factors required for polyadenylation and 3'-end cleavage of mammalian pre-mRNAs. May be responsible for the interaction of CSTF with other factors to form a stable complex on the pre-mRNA (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)