

Product datasheet for **MR206875A1V**

Mouse Crbn (BC086488) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Crbn (BC086488) AAV Particle
Tag:	Myc-DDK
Symbol:	Crbn
Synonyms:	2610203G15Rik; 2900045O07Rik; AF229032; AW108261; piL
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR206875 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGAAACCACCTGCCGCTTCTGCCTGACAGTGAAGATGAAGATGATGAAATTGAAATGGAAGTTGAAG
ACCAAGATAGTAAAGAAGCCAGAAAACCGAATATCATAAACTTTGACACCACTCTGCCAACCTCACATAC
ATACCTGGGAGCTGATATGGAGGAGTTCCACGGGAGAACTTTGCATGACGACGACAGCTGCCAGGTGATC
CCAGTCCTTCTGAGGTGCTGATGATCCTGATTCTCTGGCAGACACTCCCACTGCAGCTCTCTCACCCAC
AGGAAGTCAGCATGGTGCAGAACTTAATCCAGAAAGACAGGACCTTTGCAGTCCTTGATACAGTAATGT
GCAAGAAAGGGAAGCACAGTTTGGGACAACAGCAGAGATCTATGCCTATCGAGAAGAGCAGGAGTTTGA
ATTGAAGTAGTGAAAGTGAAGCAATTGGAAGGCAGCGTTCAAGGTCCTCGAACTTCGAACACAGTCAG
ATGGAATCCAGCAAGCTAAAGTGCAGATTTTGCCAGAGTGTGTGTTGCCGTCAACCATGTCTGCAGTGCA
GGTAGAATCACTCAATAAGTGCCAGGTATTTCTTCAAACCCATCTCTGGGAAGACCAGTATTCATGT
AAATGGTGGCAGAAATACCAGAAGAGAAAGTTTCACTGTGCAATCTAACATCATGGCCTCGTGCTGT
ATTCATTATATGATGCTGAAACATTAATGGATAGAATTAAGAAACAGCTACGTGAATGGGATGAAATCT
CAAAGATGATTCTCTTCTGAAATCCAATAGACTTTTCTTACAGAGTAGCTGCTGTCTTCTTCTATTGAT
GATGTATTGAGAATTCAGCTCCTTAAATCGGCAGTGCTATTCAACGGCTTCGCTGTGAATTGGACATCA
TGAACAAATGACTTCCCTTTGCTGTAAACAATGTCAAGAAACAGAAATAACGACAAAGAATGAAATATT
TAGTTTATCCTTATGTGGTCCAATGGCAGCATATGTGAATCCTCATGGATATGTACATGAGACACTGACT
GTGTATAAAGCGTCCAACCTGAATCTGATAGGCCGGCCTTCTACAGTGACAGCTGGTTTCCCGGTATG
CATGGACCATTGCCAGTGCAAGATCTGTGCAAGCCATATTGGATGGAATTTACAGCCACAAAAAAGA
CATGTCACCTCAAAAATTTGGGGCTTAACCTGCTCTGCTCTGTTACCCACAATTCAGAGACTGAAGAT
GAAATAAGTCCAGACAAAGTAATACTTTGTTTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR206875 protein sequence Red=Cloning site Green=Tags(s) MGNHLPLLPDSEDEDDEIEMEVEDQDSKEARKPNIINFDTSLPTSHTYLGADMEEFHGRTLHDDSCQVI PVLPEVLMILIPGQTLPLQLSHPQEVSMVRNLIQKDRTFAVLAYSNNVQEREAFQTAEIYAYREEQEFQ IEVVVKYKAIGRQRFKVLRLTQSDGIQQAQVQILPECVLPSTMSAVQVESLNKCQVFPSKPISWEDQYSC KWWQYQKRKFHCANLTSWPRWLYSLYDAETLMDRIKKQLREWDENLKDDSLPENPIDFSYRVAACLPI DVLRIQLLLKIGSAIQRLRCELDIMNKCTSLCCKQCQETETTKNEIFSLSLCGPMAAYVNPBGYVHETLT VYKASNLNLIGRPSTVHSWFGYAWTIAQCKICASHIGWKFTATKKDMSPPQKFWGLTRSALLPTIPETED EISPDKVILCL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	BC086488
ORF Size:	1293 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>BC086488</u> , <u>AAH86488</u>
RefSeq Size:	2107 bp
RefSeq ORF:	1295 bp
Locus ID:	58799
Cytogenetics:	6 E1
MW:	49.5 kDa